

| TObject | | |
|---------------------|------------------|---------------------|
| fUniqueID | kIsReferenced | kZombie |
| fBits | kHasUUID | kBitMask |
| fgDtorOnly | kCannotPick | kSingleKey |
| fgObjectStat | kNoContextMenu | kOverwrite |
| kCanDelete | kInvalidObject | kWriteDelete |
| kMustCleanup | kIsOnHeap | fgIsA |
| kObjInCanvas | kNotDeleted | |
| @~TObject | GetObjectInfo | operator new@[@] |
| MakeZombie | GetTitle | operator delete |
| DoError | HandleTimer | operator delete@[@] |
| TObject | Hash | operator delete |
| TObject | InheritsFrom | operator delete@[@] |
| operator= | InheritsFrom | SetBit |
| AppendPad | Inspect | SetBit |
| Browse | IsFolder | ResetBit |
| ClassName | IsEqual | TestBit |
| Clear | IsSortable | TestBits |
| Clone | IsOnHeap | InvertBit |
| Compare | IsZombie | Info |
| Copy | Notify | Warning |
| Delete | Is | Error |
| DistancetoPrimitive | Paint | SysError |
| Draw | Pop | Fatal |
| DrawClass | Print | AbstractMethod |
| DrawClone | Read | MayNotUse |
| Dump | RecursiveRemove | Obsolete |
| Execute | SaveAs | GetDtorOnly |
| Execute | SavePrimitive | SetDtorOnly |
| ExecuteEvent | SetDrawOption | GetObjectStat |
| FindObject | SetUniqueID | SetObjectStat |
| FindObject | UseCurrentStyle | Class |
| GetDrawOption | Write | Class_Name |
| GetUniqueID | Write | IsA |
| GetName | operator new | ShowMembers |
| GetIconName | operator new@[@] | |
| GetOption | operator new | |

| TNamed |
|--------------|
| fName |
| fTitle |
| fgIsA |
| @~TNamed |
| TNamed |
| TNamed |
| TNamed |
| operator= |
| Clear |
| Clone |
| Compare |
| Copy |
| FillBuffer |
| GetName |
| GetTitle |
| Hash |
| IsSortable |
| SetName |
| SetNameTitle |
| SetTitle |
| Is |
| Print |
| Sizeof |
| Class |
| Class_Name |
| IsA |
| ShowMembers |

| AliSpectraBothHistoManager | |
|------------------------------|----------------------------|
| fOutputList | fgIsA |
| fNRebin | |
| @~AliSpectraBothHistoManager | SetPIDHistogramByName |
| AliSpectraBothHistoManager | ManSigHistogram |
| AliSpectraBothHistoManager | ManSigHistogram |
| BookPtGenHistogram | GetNSigHistogramByName |
| BookPtGenAllChHistogram | GetVecHistogram |
| BookPtRecHistogram | GetqVecHistogram |
| BookPtRecAllChHistogram | GetVecHistogramByName |
| BookPIDHistogram | GetGenMulvsRawMulHist |
| BookNSigHistogram | GetGenMulvsRawMulHist |
| BookGenMulvsRawMulHist | GetEstgStatHist |
| BookDoubleCountsHistogram | GetOutputList |
| BookEventStatHist | SetNRebin |
| GetPtHistogram1D | GetNRebin |
| GetDCAHistogram1DMerge | |
| GetHistogram | AliSpectraBothHistoManager |
| GetHistogram1D | operator= |
| GetHistogram2D | Class |
| GetPtHistogram | Class_Name |
| GetPtHistogram | IsA |
| GetPtHistogramByName | ShowMembers |
| GetPIDHistogram | |
| GetPIDHistogram | |